

Specifications

Model	Free air displacement		Ultimate Pressure	Max. Pressure		Motor horsepower	Power Supply	Tubing needed I.D.	Weight net	Overall dimensions
	m ³ /h	l/min		mbar	Pa					
2522C-05	1.1	18.0	133	1×10 ⁶	100	93	110, 50/60	7	5.3	206×224×254
2534C-02	1.7	28.0	93	5×10 ⁵	50	93	230, 50	7	5.3	206×224×254
2546C-02	2.3	38.0	80	1×10 ⁶	100	190	230, 50	7	6.3	254×191×229
2546C-10	2.3	38.0	80	1×10 ⁶	100	190	230, 60	7	6.3	254×191×229

2522C-05, 2534C-02, 2546C-02 and 2546C-10 accessories:

- Vacuum and pressure gauge
- Vacuum and pressure regulator
- Inlet/exhaust separator jar



2522C-05



2534C-02



2546C-02 / 2546C-10

Model	Free air displacement		Ultimate Pressure	Max. Pressure		Motor horsepower	Power Supply	Tubing needed I.D.	Weight net	Overall dimensions
	m ³ /h	l/min		mbar	Pa					
2562C-02	3.4	56.7	12	—	—	250	230, 50	7	7.5	297×183×241
2562C-50										381×254×254
2567C-02	4.9	82.0	80	—	—	250	230, 50	7	7.5	297×183×241
2567C-50										381×254×254

2562C-50 and 2567C-50 accessories included with pump:

- Vacuum gauge
- Vacuum regulator
- Inlet water trap
- Exhaust silencer (option)

2562C-02 and 2567C-02 do not come with vacuum regulator, vacuum gauge and inlet water trap.



2562C-50, 2567C-50



2562C-02, 2567C-02